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

Control devices, Connection sockets, and Symbols

- A Alarm Led for the neutral plate safety circuit (red)
- B Area for memories setting (storage key  and selection/shifting keys )
- C Area of selection  and regulation  **UP/DOWN** of the Monopolar Cut/Coagulating Cut
- D Area of selection  and regulation  **UP/DOWN** of the Monopolar Coagulation
- E Area of selection  and regulation  **UP/DOWN** of the Bipolar Cut/Coagulating Cut
- F Area of selection  and regulation  **UP/DOWN** of the Bipolar Coagulation
- G Alarm Led for the HF leakage currents control circuit
- H Setting of the standard twin foot-switch pedal functioning (**PEDAL-1**) 
- I Setting of the connection sockets for the bipolar electrodes  (left socket **8a** / right socket **8b** / Both)
- 1 Socket **PEDAL-1** for for connecting the twin foot-switch pedal (standard for the activation of the monopolar or bipolar currents)
- 2 Socket **PEDAL-2** for connecting the twin foot-switch pedal (non standard, and only for the activation of the bipolar currents)
- 3 Socket **A**  for connecting monopolar hand-switch handles

- 4 Socket **B**  for connecting monopolar hand-switch handles
- 5 Socket **A**  for connecting monopolar electrodes cables with activation by foot-switch pedals
- 6 Socket **B**  for connecting monopolar electrodes cables with activation by foot-switch pedals
- 7 Socket **NEUTRAL PLATE**  for connecting the neutral electrode cable
- 8 Socket **BIPOLAR**  for connecting the bipolar electrodes cable

b1 Mains voltage selector


b2 General mains switch (green-0/I)

b3 Socket for the mains cable with fuse block

b4 Plug for the equipotential connection 

b5 Rotating regulator of the start delay of the bipolar **AUTO** coagulation (from instantaneous mode to 5sec max.)

b6 Rotating regulator of the acoustic signals intensity for the activation of the cut and coagulation

Earth protection (inside)  Alternating current

Attention: read the user manual

Apparatus of Class I Type CF, protected against the defibrillator effects (a CF type unit guarantees the highest safety level against direct and indirect contacts, notably for the allowable leakage currents). The applied part type F (floating) is protected from the earth at high and low frequencies. This kind of unit is especially indicated for direct heart application.

GEISTER MEDIZINTECHNIK GMBH guarantees safety, reliability and performances of this unit only if installation, recalibrations and repairs are carried out, using original spare parts, by personnel authorized by GEISTER and if the unit is used in compliance with given instructions in an area that meets all the applicable IEC or CEI requirements.

GEISTER is at disposal to supply, if requested, the electric diagrams and any further information needed.

This manual must be kept where the unit is employed. Please read this entire manual carefully to become familiar with each of the controls and features before making any attempt to use the equipment clinically and ask it again in case missing.

If any questions arise regarding the information contained in this manual according to specific needs, please contact the Manufacturer, directly or through the local distributor, before using the unit. In accordance with the requirements of the European Directive for medical devices 93/42 CEE and with the procedures of Company Quality System for the after-sale control of the production, the users are pleased to inform the Manufacturer about every, even little, problem of this unit.

INTRODUCTION

In a biological tissue crossed by the electric current, three effects are usually generated: Thermal, Faradic, and Electrolytic.

By using high frequency electric current, higher than 300kHz, the faradic effect is almost completely eliminated, while the electrolytic one is kept, even if with no practical result. So, the most used is the thermal one. When an electric current having such characteristics crosses with sufficient density the cellular liquid of the tissues, it warms it and generates what follows:

- 1) a heating, which is so rapid that the vapour pressure into the cells breaks their membranes and provokes their division (*pure cut*);
- 2) a heating, which is slower, and which permits to the liquid to evaporate very slowly; in this way, the coagulating parts of the tissues can coagulate (*coagulation*);
- 3) a process which is in the middle between the two phenomena described above (*coagulating cut*).

The use of the HF current presents also some risks that must be known, because their possible reduction depends also on users' behaviour:

- undesired burns on the patient tissues (i.e. where the neutral electrode is placed because of non homogeneous / insufficient contact, or in any other zone because of anomalous contacts/use of water mattresses / contacts of patient with the metallic parts of the operating table);
- undesired burns on the operator tissues (i.e. into the hand, because of an insulation leakage into the coagulation monopolar forceps);
- interferences with the functioning of other equipments (i.e. video systems) or implanted devices (pace-makers);
- **slight neuromuscular stimulations, notably with coagulation currents, both in the contact point of the active electrode and in the contact point of the neutral electrode. These stimulations are felt by patients as "electrical discharges".**

DIRECTION FOR USE

In the field of normal electro-surgery, with only the High Frequency currents, the ESU-X 300NT allows to perform all kinds of Monopolar Cut (Pure or Coagulating), Monopolar Coagulation (at low, medium, and high voltage), Bipolar Cut (Pure or Coagulating), Bipolar Coagulation (Micro, Macro, etc.), during major and medium surgery interventions into the operating theatre (open sky surgery, minimum invasive surgery, endoscopic surgery), or in any other similar place. The application medical fields of this model are the following ones:

GYNAECOLOGY, HEART SURGERY, ORTHOPAEDICS, NEUROSURGERY, OTORHINOLARYNGOLOGY, UROLOGY, MAXILLOFACIAL SURGERY, DERMATOLOGY, PLASTIC SURGERY, VASCULAR SURGERY, GENERAL AND THORACIC SURGERY, PAEDIATRIC SURGERY, EMERGENCY SURGERY, GASTROENTEROLOGY, VETERINARY, AND OTHER.

GENERAL PRECAUTIONS – *It is dangerous to ignore the following warnings!*

- 1 Every electro-surgical unit has its own characteristics and therefore, before using it, it is advisable to check its functioning, without taking into consideration the previous experiences with other devices only. Anyway, start always with very low powers, and then raise them until the required one is reached;
- 2 It is extremely dangerous to use the device if the electrical plant and the installations of the operating theatre do not comply with the current safety standards. Never use extensions for the mains cable and, if many devices are connected at the same time, ask for their compatibility to the Technical Service;
- 3 It is extremely dangerous to use accessories or instruments which are not perfectly compliant with all the applicable technical or legislative Rules, and which are not suitable for the working voltages of the device (approx. 7600Vpp "4000Vp" for the monopolar currents with crest factors equal or higher than 2; 3600Vpp "1800Vp" for the monopolar currents with crest factors lower than 2; approx. 1100Vpp "550Vp" for the bipolar currents with crest factors equal or lower than 2). Moreover, the accessories and instruments must not be old nor worn. Check always their status before the use, notably if for endoscopy. Bear in mind that:
 - All the old/worn active electrodes, accessories and cables do not work properly, and do not guarantee the perfect insulation. In addition, their unstable functioning can lead the operator to increase the output powers at dangerous levels;
 - In the user manual, for each current, the maximum output voltage "Vpp" and its variation (see the curves) according to the output power adjustment are specified. This allows the operators to choose the maximum output power that must not be overcome, in order to not exceed the rated HF insulation voltage, which is possible for each accessory;
 - The standard monopolar active electrodes for normal surgery have a stem with Ø 2.3mm (so, the standard electrode-holder handles are suitable for the electrodes having stems with this diameter).
4. Do not activate the outputs before the active electrode is in contact with the tissues, as electrical arcs can be created. They burn the tissues superficially and prevent the good effect ;
5. Keep the active electrode always clean, because otherwise it can provoke sparks or superficial carbonizations on the tissues. A dirty active electrode, an electrode in bad conditions, or an electrode with connection defects causes a reduction of the output power, as it does not have any good contact with the tissues;
6. Remember that, even if compliant with all the current standards about the electromagnetic compatibility, the unit can have interferences with other electro-medical equipments;
7. Bear in mind that, when operating on patients with pace-makers or other implanted active devices, an interference with their functioning can occur (fibrillations, etc.), and they can even be damaged (in this case, it is advisable to ask for a specific qualified advice from the Cardiology Division);
8. Never use an electro-surgical unit in presence of flammable anaesthetic gases (i.e. oxygen and nitrogen protoxide, etc.), notably when operating in cavities like thorax, abdomen, trachea, head, etc. Do not use any cleaning substances, disinfectants or flammable solvents; if used, let them evaporate before the intervention. Always remove their remaining traces from the hollow parts of the body or the cavities (i.e. umbilicus, vagina,

- etc.), and from underneath the patient. Remember that during the use a spark may cause the explosion of endogenous gases (intestine), or the fire of oxygen saturated materials (cotton, gauze, etc.);
9. Take always all the metallic objects off the patient (rings, etc.), and be sure he is not in contact with any metallic part connected to the earth, or which may conduct electricity (table, supports, etc.). Insulate with dry towels the strongly secreting parts of the body and the contacts skin-to-skin (i.e. between the arm and the body);
 10. Place always all the monitoring electrodes, which are not specifically protected, as far away as possible from the electrodes of the electrosurgical unit. It is not advisable to use needle type or very small monitoring electrodes;
 11. Use and place the neutral electrode as follows:
 - Bear in mind that, if the neutral electrodes "Split" type double section are not used, the neutral electrode safety circuit of the unit cannot control the contact between the electrode and the patient tissues (that is, it does not grant the safety of a good contact);
 - Make sure it is in perfect conditions (the worn/old neutral electrodes are extremely dangerous for the risk of the burns on the patients), and choose an area of the body as close as possible to the intervention point (the ideal would be a soft part without hairs, nor protuberant bones or superficial differences). Clean this area, shave it and massage it, in order to favour the circulation;
 - Fix it in a reliable way, without placing anything in-between, nor pressing too much, in order to avoid ischemic zones. Establish the best possible contact over the entire surface, and make sure it remains constant, especially if the patient is moved or when liquids are poured. As a matter of fact, a non homogeneous and/or insufficient contact of the neutral electrode generates both an increase of the current density in the contact points (which produces a higher temperature into the tissues, and creates burns), and a decrease of the output power into the application point (which leads the operator to raise it, dangerously);
 - Never exceed 1/3 of the max. output power for each monopolar current when using the paediatric neutral electrodes, or 1/5 when using the electrodes for babies;
 - Use the disposable neutral electrodes only once, by paying attention to the instructions on the packaging. Make sure they have the right dimensions (standard for adults with "weight body higher than 15 Kg": approx. 136cm²; standard for children with "weight body from 5 to 15 Kg": approx. 84cm²);
 - As the space between the neutral electrode and the operating area represents a sort of "path" for the HF current, be sure that it is not diagonal as regards the body, nor on the heart. Remember also that the metallic elements (prosthesis, catheters, etc.) on the path of the current may cause accumulations of current with consequent heating/burns of the surrounding tissue;
 12. Position the cables of the electrodes in a way that they do not touch the patient or any other conducting part. During the operations, place the unused active electrodes on insulating materials, far away from the patient;
 13. Always use the lowest possible power. Bear this warning in mind when intervening on patients (children or babies) for whom small neutral electrodes are used (see also point 11);
 14. Choose the bipolar technique, when operating on small portions of tissue or in cavities;
 15. Try to respect as much as possible the suggested working times, and avoid useless short-circuits between the active electrode and the neutral one;
 16. Get in contact with the Technical Service for the use of the "disposable" electrodes;
 17. When the device is switched on, check all the settings before using it on the patient, and remember that a failure can provoke an undesired increase of the power;
 18. Remember that also the use of too low powers, if combined with some particular electrodes or accessories, can cause side effects: for example, when using the Argon gas, the risk of embolism raises if the power of the spray coagulation is not able to produce quickly a rapid and impermeable eschar on the target tissues;
 19. The unit must not be used for final purposes other than those listed in this manual.

THE HF LEAKAGE CURRENTS CONTROL CIRCUIT

The unit is equipped with a leakage currents to earth control circuit because these currents represent one possible source for undesired burns on the patient or the operators (i.e. a patient who gets in touch with a metallic part of the operating table or with wet/damp towels, a patient who is placed on a water mattress for surgical needs, an operator who gets in touch with instruments or endoscopes, etc. They all are possible causes of the increase of such currents).

When the leakage currents to earth overcome 150mA (limit established by the rules), the circuit intervenes as follows:

- It automatically reduces the output power, so that the leakage currents come back within the agreed limits;
- It gives an alarm signal to the operators (red Led G-HF LEAKAGE on).

SAFETY CIRCUIT OF THE NEUTRAL ELECTRODE

The neutral electrode connection control circuit (Red Led - area A) operates in the following manners:

- 1) **With disposable/reusable electrodes with one single section (non split).** The circuit controls if the neutral electrode is connected to the cable, and if the latter is integral and correctly connected to the unit

- (socket **NEUTRAL PLATE-7**). If this is not the case, it stops the delivery of the power and gives a luminous alarm signal (the Red Led is lit, Error Code "no Np") and a buzzer (loud, intermittent);
- 2) **With disposable/reusable electrodes with twin section (split)**. The circuit works as described above at point 1, but it also checks if the quality of the contact between the electrode and the patient tissues is good enough. It operates as follows:
 - a) When the contact is optimum, the circuit does not intervene;
 - b) When the contact decreases to less than approx. 50% of the surface of a standard electrode for adults well attached, the circuit intervenes automatically reducing to max. 200W the output powers (if higher levels have been selected);
 - c) When the contact further decreases, the circuit completely stops the delivery of the power, by giving a luminous alarm signal (The Red Led is lit, Error Code "no Np") and also a buzzer (loud, intermittent).
 - 3) **When only the memories for the bipolar use (it does not require any neutral electrode) are selected (memories from 96 to 99)**, the circuit does not intervene (the 3 Led are lit, but only to show that the electrode is not connected).

INITIAL CHECKS

Initial checks

1. Make sure that the mains power supply corresponds to the technical data (see the data label on the backside of the unit), and connect the unit with the mains switch (**b2-green-** on the backside of the unit) off;
2. For a possible equipotential connection, use the plug (**b4** on the backside of the unit), and the device (**b6** on the backside of the unit) to adjust the functioning acoustic signals (the max. is clockwise). The alarm signals cannot be regulated.

DATA STORING AT THE SWITCHING ON

When switched on, or after a temporary leak of current:

The unit always stores all the settings/regulations used when switched off (see also Par. "PROGRAMS AND MEMORIES");

The unit does not keep the selection of the bipolar coagulation AUTO, and it automatically sets the coagulation **MICRO**.

The operators must select the current **AUTO** every time the device is switched on, as **the safety rules for the electro-surgical units do not allow that, when switching on, a current with automatic start/stop system is applied without having been intentionally chosen by the users.**

CONNECTION AND USE OF THE FOOT-SWITCH PEDALS

The unit is equipped with a standard twin foot-switch pedal (**34-9330**), which allows to activate or the monopolar currents or the bipolar currents (cut/coagulating cut or coagulations).

On request, the unit can also be equipped with another twin foot-switch pedal (**34-9331**), which allows the independent use of the bipolar currents, and therefore it can be very useful in the following cases:

- When surgeons want to continuously alternate the use of monopolar and bipolar currents by simply pressing on the different pedals at each time (notably in laparoscopic procedures);
- When surgeons want to use the bipolar currents through an independent pedal foot-switch, which is not the one used for the monopolar currents.

Use of the standard twin pedal foot-switch (34-9330)

Connect the pedal to the socket **PEDAL -1** (MONO/BIP) and, by pressing on the key of the area **1-H**, select the functioning mode:


- **MONO** to activate the monopolar currents of cut/coagulating cut (**yellow** pedal), or coagulation (**blue** pedal);
- **BIP** to activate the bipolar currents of cut/coagulating cut (**yellow** pedal), or coagulation (**blue** pedal), except for the coagulation **AUTO** (which is not activated by any pedal foot-switch).

Use of two twin pedal foot-switches (34-9330 and 34-9331)

Connect the pedal **34-9330** to the socket **PEDAL-1** (MONO/BIP) and the pedal **34-9331** to the socket **PEDAL-2** (BIP). When the pedal **34-9331** is connected, the Led of the area **1- H** automatically sets itself on **MONO** mode, and the unit works as follows:

- The pedal **34-9330** activates the monopolar currents of cut/coagulating cut (**yellow** pedal), or coagulation (**blue** pedal);
- The pedal **34-9331** activates the bipolar currents of cut/coagulating cut (**yellow** pedal), or coagulation (**blue** pedal), except for the coagulation **AUTO** (which is not activated by any foot-switch pedal).

FUNCTIONING, NUMBER OF USABLE ELECTRODES, AND ACTIVATION MODES

1) Monopolar functioning with HF currents

The monopolar functioning requires the use of two electrodes (an active one, and a neutral one). The current flows from the active electrode towards the neutral electrode, so that the process concerns all the tissues around the specific point where the active electrode operates.

For the use with HF currents in open sky surgery or in laparoscopy, the unit can be used:

- With 1 or 2 electrodes-holder handle/s (they can be both hand-switch type with twin key **cut/coag**, or just one can be hand-switch type with twin key **cut/coag**, while the other is normal type for use by twin pedal foot-switch **cut/coag**).

It can be used with two electrodes-holder handles (as above described) at the same time by two different operators. In order to use the unit in this way, it must be selected a coagulation mode without contact (**Fulg** or **Spray**), as stated by the International Safety Rules for the HF electro-surgical equipments (IEC 60601-2-2, par. 46.103).

When the unit is on, the bipolar modes are always possible, as specified here below at point 2).

For the use with HF currents in flexible endoscopic surgery, the unit can be used:

- With 1 flexible electrode for electro-surgery only (i.e. loop type for polypectomy) with twin pedal foot-switch **cut/coag**;

When the unit is on, the bipolar modes are always possible, as specified here below at point 2).

For the use with HF currents in endoscopic surgery under liquid in urology or gynaecology (TUR), the unit can be used:

- With 1 electrode of the monopolar resectoscope by twin pedal foot-switch **cut/coag**;

When the unit is on, the bipolar modes are always possible, as specified here below at point 2).

2) Bipolar functioning with HF currents

The bipolar functioning does not require the use of the neutral electrode, as the current flows between the tips of the bipolar electrode, and it only affects the tissues of this specific area.

The unit always allows this kind of use as follows:

- With 1 bipolar electrode (for cut/coagulating cut and/or coagulation connected to the socket **BIPOLAR** by twin pedal foot-switch **cut/coag**;
- With 1 bipolar electrode (for coagulation connected to the socket **BIPOLAR**) with automatic start/stop system (**impedance sensing**), by selecting the coagulation **AUTO**;

All the electrodes mentioned above can be used at the same time and independently, but not simultaneously, except for the two monopolar electrodes-holder handles, as specified at point 1). In this case, the first activation stops the other ones avoiding all possible errors and giving the right signal to the users (See Par. "SELF-DIAGNOSIS AND SELF-TEST" – Self-diagnosis system = Error Code uSr Act and block of the unit).

MONOPOLAR MODE FOR ELECTROSURGERY

CONNECTION AND USE OF THE ACCESSORIES

1. Socket (3-A, the led **A** above is lighted on): **hand-switch** handle (forceps for coagulation **with hand-switch for coag only**)

**** in alternative ****

Socket (5- **A**, the led **A** under is lighted on): handle (forceps for coagulation, laparoscopic instrument) for the use **by foot-switches**

(See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES")

Use only this socket to connect the instruments specified here above. If other sockets are used, a failure into the unit could occur.

Just in case of a cable with a different plug from Alsa standard type, ask for:

- **GEISTER cables**, by specifying the model and the instrument connector type
- **Adaptor (34-9590)**, for the cables with non insulated plugs Ø from 2 to 8mm, or with insulated plug Ø 4mm);

2. Socket (4- **B** the led **B** above is lighted up): second **hand-switch** handle (forceps for coagulation **with hand-switch for coag only**)

** in alternative ** (by selecting the pre-set programs 6, 7)

Socket (6- **B** the led **B** under is lighted up): second handle (forceps for coagulation) for the use **by foot-switches**

(See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES");

3. Socket (7- **NEUTRAL PLATE**): neutral electrode.

According to the connection possibilities specified here above, the operator can use:

- A Obviously, one single active electrode** (handle, forceps, instrument, etc.) **by hand-switches or foot-switches** (press on the key/pedal **yellow** for the pure cut/coagulating cut currents, or the key/pedal **blue** for the coagulation currents);
- B Alternatively, two active electrodes** (two handles, or one handle and one forceps for coagulation, or other), **both of them by hand-switches, or just one of them by hand-switches and the other by foot-switches** (press on the key/pedal **yellow** for the pure cut/coagulating cut currents, or the key/pedal **blue** for the coagulation currents)
- C At the same time, two active electrodes by two operators, both of them by hand-switches, or just one of them by hand-switches and the other by foot-switches**, (This possibility is allowed only selecting the coagulation modes without contact "**Fulg** and **Spray**", as stated by the International Safety Rules for the HF electro-surgical equipments -IEC 60601-2-2, par. 46.103);
- D Alternatively, two active electrodes** (two handles, or one handle and one forceps for coagulation, or other), **both of them by foot-switches** (This mode is possible by selecting the pre-set programs 6 and 7 "See Par. "PROGRAMS AND MEMORIES"");
- E. !!! BIPOLAR FUNCTIONING !!!**
The unit can always be used for the bipolar functioning. See Par. "BIPOLAR MODE".

In order to favour this kind of use, the unit is equipped with 3 specific pre-set programs (See Par. "PROGRAMS AND MEMORIES", programs from 3 to 5), where the following features are already memorized:

1. Activation modes and outputs available
2. Starting settings of the output powers (1 program for cardiosurgery, 1 program for neurosurgery, 1 program for general/different surgeries)

In all cases, it is possible to change all the above specified pre-settings, according to the personal need/desire, and therefore memorize these variations into the different programs (See Par. "PROGRAMS AND MEMORIES").

MONOPOLAR CURRENTS, ELECTRODES, ADJUSTMENT OF THE POWERS, ADVICE

CURRENTS FOR CUT AND COAGULATING CUT

PURE CUT (PURE) – Cut without coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key **MODE**, select the current **PURE**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the twin foot-switch, or the **yellow** key of the hand-switch handle.

Instructions for use, adjustments and electrodes

Use cutting electrodes (blade, needle, loop, Hook type, or the external part of the dissector, in laparoscopy) from 30-40W.

COAGULATING CUT (BLEND-1) – Cut with soft coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key, select the current **BLEND-1**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the foot-switch, or the **yellow** key of the hand-switch handle.

Instructions for use, adjustments and electrodes

Use cutting electrodes (blade, needle, loop, Hook type, or the external part of the dissector, in laparoscopy) from 30-40W.

COAGULATING CUT (**BLEND-2**) – Cut with a very strong spray coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key, select the current **BLEND-2**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the foot-switch, or the **yellow** key of the hand-switch handle.

Instructions for use, adjustments and electrodes

Use cutting electrodes (blade, needle, Hook type, or the external part of the dissector, in laparoscopy) from 30-40W.

COMBINED CUT (CUT ALTERNATED TO COAGULATION) (ENDO) – For flexible endoscopy with cutting phases alternated to coagulation phases

This is a constant voltage current and it is controlled by an automatic power self adjustment system (**APC System**) according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key, select the current **ENDO**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the foot-switch, or the **yellow** key of the hand-switch handle.

Instructions for use, adjustments and electrodes

Ideal only for flexible endoscopic surgery (i.e. polypectomy or papillotomy).

CURRENTS FOR COAGULATION

“FULGURATION” COAGULATION– Strong superficial sparkling effect, and optimum deep coagulating effect

This current guarantees a strong coagulating effect, both deep and superficial, and therefore it is suitable both for coagulations performed with forceps/surgical instrument, and for coagulations performed by grazing the tissues with the active electrode.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **FULG**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch , or the **blue** key of the hand-switch electrodes-holder handle. It is also possible to deliver the current by simply closing the coagulation forceps, just in case a hand-switch forceps.

Instructions for use, adjustments and electrodes

It can be used with all types of electrodes (ball, blade, or needle type, insulated coagulation forceps, loop electrodes in conization,

Hook electrodes, dissectors, and forceps in laparoscopy) from 40-50W.

“PIN POINT” COAGULATION – Medium-Low superficial sparkling effect, and optimum deep coagulating effect

This current guarantees a strong deep coagulating effect and a normal superficial coagulating effect. It is suitable both for coagulations performed with forceps/surgical instrument and for coagulations performed directly with the active electrode, just in case the operators prefer a superficial effect which is less strong than the one of the FULG coagulation.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **PIN POINT**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch, or the **blue** key of the hand-switch electrodes-holder handle. It is also possible to deliver the current by simply closing the coagulation forceps, just in case a hand-switch forceps is used.

Instructions for use, adjustments and electrodes

It can be used with all types of electrodes (ball, blade, or needle type, insulated coagulation forceps, loop electrodes in conization,

Hook electrodes, dissectors, and forceps in laparoscopy) from 40-50W.

“SOFT” COAGULATION – Low superficial sparkling effect, and good deep coagulating effect

This current guarantees a good deep coagulating effect, and a very low superficial coagulating effect.. It is suitable for coagulations performed or by forceps/surgical instrument or touching the tissues by an active ball electrode of 4-5mm at least.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **SOFT**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch, or the **blue** key of the hand-switch electrodes-holder handle. It is also possible to deliver the current by simply closing the coagulation forceps, just in case a hand-switch forceps is used.

Instructions for use, adjustments and electrodes

Use coagulation electrodes (ball type, coagulation insulated forceps in laparoscopy) from 50-60W.

“SPRAY” COAGULATION – Very strong superficial sparkling effect, and good deep coagulating effect

This current guarantees a good deep coagulating effect, and a very strong superficial coagulating effect. It is suitable both for the coagulations performed with forceps/surgical instrument, and for coagulations performed directly with the active electrode, even far away from the tissues.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **SPRAY**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch, or the **blue** key of the hand-switch electrodes-holder handle. It is also possible to deliver the current by simply closing the coagulation forceps, just in case a hand-switch forceps is used.

Instructions for use, adjustments and electrodes

Use cutting and coagulation electrodes (blade, needle, ball type, and coagulation insulated forceps, in laparoscopy,

conization electrodes in gynaecology) from 40-50W.

ADVICE

OPEN SKY SURGERY

For the Cut, use as follows:

Current **PURE**, from 40-50W, for a cut without the coagulating effect (it is not very common in surgery, where the coagulating cut is preferred);

Current **BLEND-1**, from 40-50W, for a cut with a medium coagulating effect;

Current **BLEND-2**, from 40-50W, for a cut with a very strong superficial coagulating effect, *spray* type;

For the Coagulation, use as follows:

Current **FULG**, from 40-50W, for superficial coagulations or deep coagulations with the active electrode or the forceps. It is the most suitable coagulation mode when a forceps or an electrode is used;

Current **PIN POINT**, from 40-50W, for superficial coagulations or deep coagulations with the active electrode or the forceps, if a softer superficial effect is desired, compared to the **FULG** coagulation;

Current **SOFT**, from 40-50W, for superficial or deep coagulations, which are softer than those of the **PIN POINT** coagulation, which is not suitable at all for the use with small electrodes with which it tends to cut;

Current **SPRAY**, from 40-50W, for a very strong superficial coagulation with the active electrode (without any contact with the tissues), and a good deep coagulation with the forceps (it is the ideal coagulation mode for the operators who prefer the active electrode to coagulate).

LAPAROSCOPIC SURGERY

For the Cut, use as follows:

See above "OPEN SKY SURGERY" **PURE**, **BLEND-1**, **BLEND-2**. The current **BLEND-2** is very efficacious when using the Hooks or the external part of the forceps/dissectors to cut with a strong coagulating effect.

For the Coagulation, use as follows:

See above "OPEN SKY SURGERY" **FULG**, **PIN POINT**, **SPRAY**. The current **SPRAY** is very efficacious when using the Hooks or the external part of the forceps/dissectors to coagulate.

FUNCTIONING FOR THE FLEXIBLE ENDOSCOPIC SURGERY

In order to favour this kind of use, the unit is equipped with 1 specific pre-set program (See Par. "**PROGRAMS AND MEMORIES**", program 8), where the following features are already memorized:

1. Activation mode and output available
2. In this program are not set output powers or currents for the normal monopolar use because this choice depends on the operator.
3. In all cases, it is possible to set standard pre-set currents/output powers and therefore memorize these settings into the program (See Par. "**PROGRAMS AND MEMORIES**").

CONNECTION AND USE OF THE ACCESSORIES

1. Socket (5-A, the led **A** under is lighted on): flexible electrode for electro-surgery with pedal foot-switch (See Par. "**CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES**").

Use only this socket to connect the above mentioned instruments. The use of different sockets provokes a damage to the unit.

Just in case, ask for:

The adaptors for the cables with different plug:

Adaptor- 34-9590, for cables with non insulated plugs Ø from 2 to 8mm, or with insulated plug Ø 4mm;

2. Socket (7- **NEUTRAL PLATE**): neutral electrode.

According to the connection possibilities listed above, the unit allows to:

- A. Use the flexible electrode for the normal HF electro-surgery by foot-switch activation (See Par. "**CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES**"). Press on the **yellow** pedal for the pure/coagulating cut currents, and the **blue** pedal for the coagulation currents.
- C. **!!! BIPOLAR FUNCTIONING!!!**
The bipolar functioning is always possible. See Par. "**BIPOLAR MODE**".

MONOPOLAR CURRENTS, ELECTRODES, ADJUSTMENT OF THE POWERS, ADVICE

CURRENTS FOR CUT

PURE CUT – Cut without coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key **MODE**, select the current **PURE**, then adjust the power with the regulation keys **UP/DOWN**;
- **For the delivery (loud acoustic signal, and yellow Led), press on the yellow pedal of the twin foot-switch**

Instructions for use, adjustments and electrodes

Use the flexible electrodes of the endoscopes, from 20-30W for the papillotomy, from 20-30W for the polypectomy with polypus

Ø 5mm, and from 40-50W for polypus Ø 6mm or more.

COAGULATING CUT (BLEND-1) – Cut with soft coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key, select the current **BLEND-1**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the twin foot-switch.

Instructions for use, adjustments and electrodes

Use the flexible electrodes of the endoscopes, from 20-30W for the papillotomy, from 20-30W for the polypectomy with polypus

Ø 5mm, and from 40-50W for polypus Ø 6mm or more.

COAGULATING CUT (BLEND-2) – Not suitable at all.

COMBINED CUT (CUT ALTERNATED TO COAGULATION) (ENDO) – For flexible endoscopy with cutting phases alternated to coagulation phases

This is a constant voltage current and it is controlled by an automatic power self adjustment system (APC System) according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key, select the current **ENDO**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the twin foot-switch.

Instructions for use, adjustments and electrodes

Use the flexible electrodes of the endoscopes, from 70-80W for the papillotomy, from 70-80W for the polypectomy with polypus

Ø 5mm, and from 90-100W for polypus Ø 6mm or more.

CURRENTS FOR COAGULATION

“FULGURATION” COAGULATION AT HIGH VOLTAGE - Strong superficial sparkling effect, and optimum deep coagulating effect

This current guarantees a strong coagulating effect, both deep and superficial, and therefore it is suitable both for coagulations performed with forceps/surgical instrument, and for coagulations performed by grazing the tissues with the active electrode.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **FULG**, then adjust the power with the regulation keys **UP/DOWN**;
- *For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch.*

Instructions for use, adjustments and electrodes

Use the flexible electrodes of the endoscopes, from 20-30W.

“PIN POINT” COAGULATION AT MEDIUM VOLTAGE – Medium-Low superficial sparkling effect, and optimum deep coagulating effect

This current guarantees a strong deep coagulating effect, and a normal superficial coagulating effect. It is suitable both for coagulations performed with forceps/surgical instrument, and for coagulations performed by grazing the tissues with the active electrode, just in case the operators prefer a superficial effect which is less strong than the one of the FULG coagulation.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **PIN POINT**, then adjust the power with the regulation keys **UP/DOWN**;
- *For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch.*

Instructions for use, adjustments and electrodes

Use the flexible electrodes of the endoscopes, from 20-30W.

“SOFT” COAGULATION AT LOW VOLTAGE – Low superficial sparkling effect, and good deep coagulating effect

This current guarantees a good deep coagulating effect, and a very low superficial coagulating effect. It is suitable for coagulations performed or by forceps/surgical instrument or by touching the tissues with an active ball electrode of 4-5mm at least.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **SOFT**, then adjust the power with the regulation keys **UP/DOWN**;
- *For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch.*

Instructions for use, adjustments and electrodes

Use the flexible electrodes of the endoscopes, from 30-40W.

“SPRAY” COAGULATION AT VERY HIGH VOLTAGE – Very strong superficial sparkling effect, and good deep coagulating effect

This current guarantees a good deep coagulating effect, and a very strong superficial coagulating effect. It is suitable both for the coagulations performed with forceps/surgical instrument, and for coagulations performed directly with the active electrode, even far away from the tissues.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **SPRAY**, then adjust the power with the regulation keys **UP/DOWN**;
- *For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch.*

Instructions for use, adjustments and electrodes

Use the flexible electrodes of the endoscopes, from 20-30W.

ADVICE

For the Cut, use as follows:

Current **PURE**, for a cut without coagulating effect and constant delivery of the power (from 20-30W for papillotomy or polypectomy with polypus Ø until 5mm, and from 30-40W with polypus Ø 6mm or more);

Current **BLEND-1**, for a cut with a soft coagulating effect and constant delivery of the power (from 20-30W for papillotomy or polypectomy with polypus Ø until 5mm, and from 30-40W with polypus Ø 6mm or more);

Current **ENDO**, for a cut-coagulation combination (alternated phases of cut and coagulation), and self regulation of the power (from 40-50W for papillotomy or polypectomy with polypus Ø until 5mm, and from 70-80W with polypus Ø 6mm or more).

For the Coagulation, use as follows:

Current **FULG**, from 30-40W;

Current **PIN POINT**, from 40-50W, for a superficial effect which is softer than the one of the **FULG** coagulation.

FUNCTIONING FOR THE UNDER LIQUID ENDOSCOPIC SURGERY (TUR, TURP, TURV, etc.)

In order to favour this kind of use, the unit is equipped with 1 specific pre-set program (See Par. "PROGRAMS AND MEMORIES", programs 8), where the following features are already memorized:

1. Activation mode and output available for the connection of the resectoscope
 2. In this program are not set output powers or currents because this choice depends on the operator.
- In all cases, it is possible to set standard pre-set currents/output powers and therefore memorize these settings into the different programs (See Par. "PROGRAMS AND MEMORIES").**

CONNECTION AND USE OF THE ACCESSORIES

A) Socket (5-A, the led A under is lighted on): monopolar rigid resectoscope with pedal foot-switch (See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES").

Use only this socket to connect the above mentioned instruments. The use of different sockets provokes a damage to the unit.

Just in case, ask for:

The adaptors for the cables with different plug:

Adaptor- 34-9590, for cables with non insulated plugs Ø from 2 to 8mm, or with insulated plug Ø 4mm;

B) Socket (7- NEUTRAL PLATE): neutral electrode.

According to the connection possibilities listed above, the unit allows to:

A. Use the resectoscope by foot-switch activation (See PAR. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES").

Just press on the **yellow** pedal for the pure cutting/coagulating cutting currents, and on the **blue** pedal for the coagulation currents.

B. !!! BIPOLAR FUNCTIONING!!!

The bipolar functioning is always possible. See Par. "BIPOLAR MODE".

MONOPOLAR CURRENTS, ELECTRODES, ADJUSTMENT OF THE POWERS, ADVICE

CURRENTS FOR CUT

PURE CUT (PURE) – Cut without coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key **MODE**, select the current **PURE**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the twin foot-switch.

Instructions for use, adjustments, and electrodes

Use the specific electrodes for the resectoscope from:

70-80W with wire loop electrodes, for soft tissues (i.e. into the bladder);

90-100W with wire loop electrodes, for strong tissues (i.e. TUR);

160-180W with ribbon loop (cutting/vaporization) electrodes, for soft tissues (i.e. into the bladder);

200-210W with ribbon loop (cutting/vaporization) electrodes, for strong tissues (i.e. TUR);

100-120W with small rolls for vaporisation;

200-210W with big rolls for vaporisation.

COAGULATING CUT (BLEND-1) – Cut with soft coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area C):

- With the selection key **MODE**, select the current **BLEND-1**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the twin foot-switch pedal

**Instructions for use, adjustments, and electrodes
As it is detailed for the current *PURE*.**

COAGULATING CUT (BLEND-2) – Not suitable at all.

COMBINATION CUT (CUT ALTERNATED TO COAGULATION) (ENDO) – Not suitable at all.

CURRENTS FOR COAGULATION

“FULGURATION” COAGULATION AT HIGH VOLTAGE - Strong superficial sparkling effect, and optimum deep coagulating effect

This current guarantees a strong coagulating effect, both deep and superficial. It is suitable both for coagulations performed with forceps/surgical instrument, and for coagulations performed by grazing the tissues with the active electrode.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key, select the current **FULG**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch.

Instructions for use, adjustments, and electrodes

Use the loop/ball electrodes for the resectoscope from 60-70W.

“PIN POINT” COAGULATION AT MEDIUM VOLTAGE – Medium-Low superficial sparkling effect, and optimum deep coagulating effect

This current guarantees a strong deep coagulating effect and a normal superficial coagulating effect. It is suitable both for coagulations performed with forceps/surgical instrument, and for coagulations performed by grazing the tissues with the active electrode, just in case the operators prefer a superficial effect which is less strong than the one of the FULG coagulation.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key **MODE**, select the current **PIN POINT**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch .

Instructions for use, adjustments, and electrodes

Use the loop/ball electrodes for the resectoscope from 60-70W.

“SOFT” COAGULATION AT LOW VOLTAGE – Not suitable.

“SPRAY” COAGULATION AT VERY HIGH COLTAGE (*SPRAY*) – Very strong superficial sparkling effect, and good deep coagulating effect

This current guarantees a good deep coagulating effect, and a very strong superficial coagulating effect. It is suitable both for the coagulations performed with forceps/surgical instrument, and for coagulations performed directly with the active electrode, even far away from the tissues.

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area D):

- With the selection key **MODE**, select the current **SPRAY**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch.

Instructions for use, adjustments, and electrodes

Use the loop/ball electrodes for the resectoscope from 60-70W.

ADVICE

For the Cut with wire loop electrodes, use as follows:

Current **PURE**, from 70-80W into the bladder or on soft tissues, for a cut without coagulating effect;
Current **BLEND-1**, from 70-80W into the bladder or on soft tissues, for a cut with soft coagulating effect;
Current **PURE**, from 90-100W on the prostate, for a cut without coagulating effect;
Current **BLEND-1**, from 90-100W on the prostate, for a cut with soft coagulating effect.

For the Cut with ribbon loop electrodes, use as follows:

Current **PURE**, from 160-180W into the bladder or on soft tissues, for a coagulating cut;
Current **BLEND-1**, from 160-180W into the bladder or on soft tissues, for a cut with stronger coagulating effect;
Current **PURE**, from 200-210W on the prostate, for a coagulating cut;
Current **BLEND-1**, from 200-210W on the prostate, for a cut with stronger coagulating effect.

NB: In order to optimize the functioning and obtain the perfect cut without any “sticking” of the tissues on the electrode, the above mentioned powers maybe have to be changed through small variations of 10-15W each time.

For the Vaporisation with rolls, use as follows:

Current **PURE**, from 100-120W with small rolls, and from 200-210W with big rolls.

For the Coagulation, use as follows:

Current **SPRAY** from 40-50W. Optimum to coagulate or by loop/blade/ball electrodes, even without any contact with the tissues;

Current **FULG**, from 60-70W, for a coagulation by loop/blade/ball electrodes.

Current **PIN POINT**, from 60-70W, for a coagulation by loop/blade/ball electrodes with contact between the electrode and the tissues.

BIPOLAR MODE

INTRODUCTION

As far as the bipolar functions, please note as follows:

- They are always available, as alternative to all the other monopolar functioning modes;
- They can be used alone, without connecting the neutral electrode. See. the pre-set programs in next paragraph;
- Besides the standard twin foot-switch (**93-9330**, for the monopolar or bipolar use), the unit can also be equipped or with another twin foot-switch (**34-9331**, for the bipolar use only), on request. (See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES");
- The unit is equipped with a coagulating current **AUTO** (impedance sensing type) with automatic start/stop system (switching on with adjustable delay from 0.5 to 5sec by the control device **b5** which is placed on the backside of the unit, when the tips of the forceps get in touch with the non coagulated tissues – switching off when the tissues are coagulated). This system does not require any pedal foot-switch;

CONNECTION AND USE OF THE ACCESSORIES

The unit is equipped with one socket for the connection of the bipolar electrodes (**8-BIPOLAR**).

1. Socket **8- BIPOLAR-A** : bipolar electrode (connected cable without any particular polarity) for use **by foot-switch** for cutting and/or coagulation. The foot-switch for the **AUTO** coagulation with automatic start/stop system is not used;

In order to facilitate the bipolar use, the unit is equipped with 1 pre-set program (See Par. "PROGRAMS AND MEMORIES",

program 9) where the following features are already memorized :

- Disconnection of the neutral electrode safety circuit, which is not required for the bipolar use only (as a matter of fact, through the 3 lit alarm Led, the circuit just informs the operators that the neutral electrode is not connected);
- Impossibility to use the monopolar currents (the displays for the monopolar currents powers are switched off);
- Activation modes and outputs available according to the different settings

Just in case of cables with connection plugs which are different from GEISTER standard type, ask for:

- **GEISTER cables**, by specifying the model and the instrument connector type;
- **Adaptor, 34-9591V** for cables with connector (Valleylab type).
- **Adaptor, 34-9591E** for cables with connector (coaxial type, Erbe or Storz type).
- **Adaptor, 34-9591M** for cables with connector (coaxial type, Martin or Bertchold type).

BIPOLAR CURRENTS, ELECTRODES, ADJUSTMENT OF THE POWERS, ADVICE

CURRENTS FOR CUT

PURE CUT– Cut with minimum coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area L):

- With the selection key **MODE**, select the current **PURE**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the twin foot-switch.

Instructions for use, adjustments, and electrodes

Use the specific electrodes, from:

40-50W, with cutting electrodes, needle type;

COAGULATING CUT – Cut with coagulating effect

This current is controlled by an automatic system (ADC System) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues.

Control devices to be used (area L):

- With the selection key **MODE**, select the current **BLEND**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**loud** acoustic signal, and **yellow** Led), press on the **yellow** pedal of the twin foot-switch.

Instructions for use, adjustments, and electrodes

As it is detailed for current PURE.

CURRENTS FOR COAGULATION

“MICRO” COAGULATION – Very soft coagulating effect

This is a constant voltage current and it is controlled by an automatic power self adjustment system (**APC System**) according to the characteristics of the tissues. It is the most suitable current when operating on extremely delicate tissues (i.e for neurosurgery).

Control devices to be used (area F):

- With the selection key **MODE**, select the current **MICRO**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch.

Instructions for use, adjustments, and electrodes

Use coagulation forceps with tips of: 0.5mm (from 0.5-1W), 1mm (from 1-2W), and 2mm (from 2-4W);

Use coagulation forceps for laparoscopy, from 30-40W.

“AUTO” COAGULATION – Very soft coagulating effect and automatic start/stop system (impedance sensing type)

This current is the same as the MICRO current, as far as the performances and the indications are concerned.

Control devices to be used (area F):

- With the selection key **MODE**, select the current **AUTO**, then adjust the power with the regulation key **UP/DOWN** s;
- For the delivery (**acute** acoustic signal, and **blue** Led), touch the tissues with the forceps tips after having regulated the starting delay (from 0.5 to 5sec) by turning the control device **b5** on the backside of the unit (0.5sec on the left, 5sec on the right);
- For the switching off it is possible either to disconnect the forceps from the tissues or wait until the current disconnects automatically, once the tissues are coagulated

For safety reasons, the unit is also provided with an automatic stop control, which is activated after approximately 6 seconds of functioning (if the forceps remains on the tissues, it starts again after a stop interval which is exactly the same as the delay set for the start).

Instructions for use, adjustments, and electrodes

Use coagulation forceps with tips of: 0.5mm (from 0.5-1W), 1mm (from 1-2W), and 2mm (from 2-4W);

Not suitable for the use in laparoscopy.

“MACRO” COAGULATION – Strong coagulating effect

This current is controlled by an automatic system (**ADC System**) for the self regulation of the functioning, which guarantees that the power selected by the operators remains constant according to the characteristics of the tissues. It is the most suitable current for a strong coagulating effect on big dimensions electrodes, and therefore it is ideal both for the operators who desire a bipolar coagulation with an effect which is “similar” to the monopolar coagulation, and for those who want to coagulate in laparoscopy or on big tissues areas.

Control devices to be used (area F):

- With the selection key **MODE**, select the current **MACRO**, then adjust the power with the regulation keys **UP/DOWN**;
- For the delivery (**acute** acoustic signal, and **blue** Led), press on the **blue** pedal of the twin foot-switch .

Instructions for use, adjustments, and electrodes

Use coagulation forceps with tips of: 1mm (from 5-10W), 2mm (from 7-15W);

Use coagulation forceps for laparoscopy, from 40-50W;

ADVICE

For the Cut:

In order to optimize the functioning of the bipolar cut and obtain the perfect cut without any “sticking” of the tissues on the electrode, the above mentioned powers maybe have to be changed through small variations of 10/15W each time.

For the Coagulation:

When bipolar forceps are used, it is advisable to keep their tips damp, for example with physiologic solution, in order to reduce the omnipresent phenomenon of the “sticking” of the tissues. It is also possible to use forceps with “anti-sticking” tips, or better with irrigation system.

PROGRAMS AND MEMORIES

The unit is provided with 10 memories (9 programmable memories and 1 for the storage of the data at the switching off), which allow to set 9 complete functioning programs, that can be selected by simply pressing the shifting keys UP/DOWN until the desired memory is chosen. 2 memories are totally free, and can be set by the operators, while the remaining ones are pre-set in order to facilitate the use of the unit:

How it must be set a free memory

1. Select the desired memory with the selection keys UP/DOWN of the area **B-MEMORY**. (It is possible to shift the memories from 1, 2, 3 and so on until 9 with the key UP, but also to shift them in the opposite direction, from 9 to 1, with the key DOWN);
2. Make (in the area **H**) all the adjustments related to foot-switch activation (by the standard twin foot-switch **34-9330** for monopolar/ bipolar use or by the independent twin foot-switch **34-9331** for the bipolar use only), except for the **AUTO** bipolar coagulation;
3. Choose all the desired current/s and the adjustment of their power/s.
For example: as the unit is equipped with 4 different monopolar coagulation currents, it is possible to set a memory, through the keys of the area **D**, with one single current and its power, or with four currents and their different powers which can be then changed during the use, by simply pressing the selection key.
In the first case only the desired current must be selected by the selection key and only its power must be set by the adjustment keys.
In the second case, one upon the other, each current and their power must be set by the selection key and by the adjustment keys;
4. After having set all the parameters, the program has to be confirmed by pressing the storage key **MEM** of the area **B**.

Is it possible to modify during the use the adjustments of the powers already set in a memory?

Yes, it is always possible, through the adjustment keys of the power of the used current, and this variation is shown by the number of the memory which blinks. The number stops blinking if:

- The storage key **MEM** of the area **B** is pressed, in order to confirm the variation done;
- One of the two selection keys **UP/DOWN** of the area **B** is pressed, in order to come back to the powers and/or the functioning mode which was set before the variation.

Is it possible to modify the adjustments of the power/s and/or the selected settings of a memory before the use?

Yes, it is always possible, through the different adjustment/selection keys, and this variation is shown by the number of the memory which blinks. The number stops blinking if:

- The storage key **MEM** of the area **B** is pressed, in order to confirm the variation done;
- One of the two selection keys **UP/DOWN** of the area **B** is pressed, in order to come back to the powers and/or the functioning mode which was set before the variation.

Is it possible to come back to the original power and/or functioning settings, if changed but not confirmed (while the number of the memory is blinking)?

Yes, it is. One of the two selection keys **UP/DOWN** of the area **B** must be pressed (the number of the memory stops blinking, and the original data are set again).

Is it possible to use the unit without setting any memory?

Yes, it is, but it is less practical, as each type of intervention may require different parameters.

What is the best way to use a preset program?

A preset program helps the users assuring a specific setting needed or required for performing a specific surgical procedure (the setting changes according to the procedure or according to the behaviour of each operator) and, because of this reason, it is possible to set and memorize a starting program. Then, if during the use, it is needed a changing (Current or Power can be obviously changed as above specified), the best way is:

- The following the operation using the different selection (the number of memory blinks).
- The resetting the program, at the end of the operation, to the original data by pushing, without distinction, the selection keys **UP/DOWN** of the area **B**

What happens if the unit is switched off while the memory number is blinking?

At the moment of the switching on the unit starts with the setting used when switched off. In this case, it means with the same memory (The number still blinking) and with the same setting.

PRE-SET MEMORIES (nr. 3) – Use in neurosurgery

The pre-setting concerns:

1. Activation modes and outputs available.
2. Starting settings of the output powers for the normal HF neurosurgery when the Argon section switched off) (70W for all the monopolar cutting currents **PURE**, **BLEND 1**, **BLEND 2**, 40W for the monopolar coagulating currents **FORCED**, **PIN POINT**, **SPRAY**, 20W for the more delicate **MICRO** bipolar coagulating current)

The activation modes and the outputs available can not be changed, but all the currents and the output powers can obviously be changed by setting the desired levels and by pressing the memorisation key MEM of the area B.

The pre-setting of these memories allows to use the unit as follows:

A) With 2 active electrodes:

One active electrode or by hand-switches (connected to the socket **3-A**) or by foot-switches (connected to the socket **5-A**).

In order to activate the different active electrodes shall be pressed (on the handles or on the twin foot switch) the **yellow** key/pedal for the pure/coagulating cut and the **blue** key/pedal for the coagulation currents.

In order to use two electrodes-holder handles (as above described) at the same time by two different operators, it must be selected a coagulation mode without contact (**Fulg** or **Spray**), as stated by the International Safety Rules for the HF electro-surgical equipments (IEC 60601-2-2, par. 46.103).

Another active electrode, for the normal electro-surgery by hand-switches (connected to the socket **4-B**) pressing the **yellow** key for the pure/coagulating cut and the **blue** key for the coagulation currents (See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES");

C. !!! BIPOLAR FUNCTIONING !!!

The bipolar functioning is always available. See Par. "BIPOLAR MODE".

PRE-SET MEMORIES (nr.4) – Use in cardio surgery.

The pre-setting concerns:

1. Activation modes and outputs available.
2. Starting settings of the output powers for cardio-surgery (70W for all the monopolar cutting currents **PURE**, **BLEND 1**, **BLEND 2**, 70W for the monopolar coagulating currents **FORCED**, **PIN POINT**, **SPRAY** setting itself on **SPRAY** because it is the more appreciated by the cardio surgeons)
3. Activation modes and outputs available for cardio-surgery.
4. Starting settings when the Argon section is switched on (Automatic selection of the **SPRAY** coagulation current that is needed for the Argon coagulation "the Led in the area **D** blinks when this type of current is not the one used also for the electrosurgical mode", automatic delivery of the **SPRAY** current at 60W when the **blue** key on the electrodes-holder handle for Argon gas (**34-9730**) is pressed, automatic setting and delivery of the Argon gas flow at the level of "5-6 lt/min").

The activation modes and the outputs available can not be changed, but all the currents, the output powers and the Argon gas flow (if it is used) can obviously be changed by setting the desired levels and by pressing the memorisation key MEM of the area B.

The pre-setting of these memories allows to use the unit as follows:

A) With 2 active electrodes:

One active electrode or by hand-switches (connected to the socket **3-A**) or by foot-switches (connected to the socket **5-A**).

In order to activate the different active electrodes shall be pressed (on the handles or on the twin foot switch) the **yellow** key/pedal for the pure/coagulating cut and the **blue** key/pedal for the coagulation currents.

In order to use two electrodes-holder handles (as above described) at the same time by two different operators, it must be selected a coagulation mode without contact (**Fulg** or **Spray**), as stated by the International Safety Rules for the HF electro-surgical equipments (IEC 60601-2-2, par. 46.103).

Another active electrode, for the normal electro-surgery by hand-switches (connected to the socket **4-B**) pressing the **yellow** key for the pure/coagulating cut and the **blue** key for the coagulation currents.(See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES");

C. !!! BIPOLAR FUNCTIONING !!!

The bipolar functioning is always available. See Par. "BIPOLAR MODE".

PRE-SET MEMORIES (nr. 5) – Use in open sky/laparoscopic surgery.

The pre-setting concerns:

1. Activation modes and outputs available.
2. Starting settings of the output powers (80W for all the monopolar cutting currents and 70W for all the monopolar coagulating currents)
3. Activation modes and outputs available.

The activation modes and the outputs available can not be changed, but all the currents and the output powers can obviously be changed by setting the desired levels and by pressing the memorisation key MEM of the area B.

The pre-setting of these memories allows to use the unit as follows:

A) With 2 active electrodes:

One active electrode or by hand-switches (connected to the socket **3-A**) or by foot-switches (connected to the socket **5-A**).

In order to activate the different active electrodes shall be pressed (on the handles or on the twin foot switch) the **yellow** key/pedal for the pure/coagulating cut and the **blue** key/pedal for the coagulation currents.

In order to use two electrodes-holder handles (as above described) at the same time by two different operators, it must be selected a coagulation mode without contact (**Fulg** or **Spray**), as stated by the International Safety Rules for the HF electro-surgical equipments (IEC 60601-2-2, par. 46.103).

Another active electrode, for the normal electro-surgery by hand-switches (connected to the socket **4-B**) pressing the **yellow** key for the pure/coagulating cut and the **blue** key for the coagulation currents.(See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES");

C. !!! BIPOLAR FUNCTIONING !!!

The bipolar functioning is always available. See Par. "BIPOLAR MODE".

PRE-SET MEMORIES (nr.6 and 7) – These two programs allow to use, always by the same twin foot-switch, two different monopolar active electrodes (one blade or needle electrode for cutting inserted in a handle and one forceps for coagulation with connection cable to the unit) simply setting the same currents or output powers and alternating the two programs (i.e. the program 89 or for using the cutting or coagulating current "instrument for cutting or coagulating cutting" and the program 90 for using the coagulating current "instrument for coagulation")

The pre-setting does not concern the currents or their powers which have to be set according to the operators needs and memorised as specified in the first part of the paragraph (for this kind of use, it is advisable to select and set the same currents and the same powers in both memories).

The pre-setting of these memories allows to use the unit as follows:

- A. **Selecting the memory 89, one electrodes-holder handle** (connected to the socket **5-A**) can be used by the foot-switch pedal. (See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES") pressing (on the twin foot switch) the **yellow** pedal for the pure/coagulating cut and the **blue** pedal for the coagulation currents;
- B. **Selecting the memory 90, one coagulation forceps** (connected to the socket **6**) can be used by the foot-switch pedal. (See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES") pressing (on the twin foot switch) the **blue** pedal for the coagulation currents;

The activation modes and the outputs available can not be changed, but it is possible to set pre-set currents/output powers and memorize the settings into these programs by pressing the key MEM of the area B (See the first part of this paragraph).

PRE-SET MEMORIES (nr.8) – Use in flexible endoscopic surgery or in endoscopic surgery under liquid.

The pre-setting concerns:

1. Activation mode and output available.
2. In these programs are not set output powers or currents because this choice depends on the operator.

The activation modes and the outputs available can not be changed, but it is possible to set standard currents/output powers and memorize the settings into these programs by pressing the key MEM of the area B (See the first part of this paragraph).

The pre-setting of these memories allows to use the unit as follows:

- A. With the flexible endoscope or with the rigid resectoscope in (connected to the socket **5-A**) by pedal foot-switches (See Par. "CONNECTION AND USE OF THE PEDAL FOOT-SWITCHES"), pressing (on the twin foot switch) the **yellow** pedal for the pure/coagulating cut and the **blue** pedal for the coagulation currents;
- B. **!!! BIPOLAR FUNCTIONING !!!**

The bipolar functioning is always available even if it is not used in the endoscopy under liquid but only, hardly ever, in flexible endoscopy. See Par. "BIPOLAR MODE".

PRE-SET MEMORIES (nr. 9) – Use in bipolar mode only.

The pre-setting does not concern the current types that have to be used, nor their powers, which have to be set according to the operators' needs, and therefore memorised as shown at points 3 and 4 of the previous paragraph.

The pre-setting of these memories allows to use the unit as follows:

- A. Without connecting the neutral electrode, as it is not required (through the alarm Led, the safety circuit just informs the operators that the neutral electrode is not connected);
- B. Without using the monopolar currents (their displays are switched off);

It is possible to set standard programs including: the currents and their powers, the number of twin foot-switches (1 or 2) and memorize the settings into the programs by pressing the key MEM of the area B (See the first part of this paragraph).

SELF-DIAGNOSIS and SELF TEST SYSTEM, CONTROL of the ERROR CODES

The unit is equipped with a self-diagnosis system, which automatically intervenes in case of functioning problems, included the delivery of powers if are higher than the set ones. It stops the device and informs the operators through specific alarm signals (*Error Codes*), which can be either acoustic (loud and intermittent sound) or also visual (see the Tables *Codes for Wrong Use or Activation Circuit Fault – Codes for System Fault*).

When the device is turned on, this control system carries out a complete self-test cycle over the hardware and the software (all the Led and displays light up, and the loudspeaker is activated). If everything works properly, this phase ends by a short acoustic signal.

If the system detects any problem, the faults are indicated as described above, and the operator must either try to eliminate all the possible causes, or switch the equipment off and ask for a technical assistance.

When the device is running, this control system checks the functioning of the unit, included the output powers (2000 times per second), by repeating the complete test on the functioning every 20 minutes (in an activation interval).

Even in this case, the problems are indicated as described above ("*When the device is turned on*").

Control of the Error Codes

The unit stores the last 32 Error Codes that have been detected, in order to help the operators in finding a solution for the problems/faults. The procedure to follow for this control is detailed into the Service Manual.

Inside the E²PROM hexadecimal value corresponding to one of error codes reported on "Error code" column in the above table are stored.

Error description	Error code	Stored code in E ² PROM (select "0d" in CAL mode) ¹
Codes for Wrong Use or Activation Circuits Fault		
For all the "U" problems, the operators can act by eliminating the cause of the intervention of the self-diagnosis system. In all the other cases, ask for the technical assistance, after having checked the signal by switching off and on the unit		
U) Problem into the neutral electrode (broken cable, non connected cable to the unit or to the neutral electrode, neutral electrode "split" type not attached to the tissues)	no nP	-
Failure into the neutral electrode control circuit (internal fault)	Err nPC	81
Failure into the control circuit of the contact between the neutral electrode/the patient (internal fault)	Err 53	35
Failure of an activation hand-switch at the switching on	Err Hnd	82
U) Activation hand-switch pressed at the switching on	Err PEd	83
Failure of a pedal foot-switch at the switching on	Err PEd	83
U) Pedal foot-switch pressed at the switching on	USr Act	-
U) Contemporary activation not allowed of the two activation devices	USr Act	-
U) Key on the control panel pressed at the switching on	Err 14	0E
Failure of a key on the control panel	Err 14	0E
Codes for System Fault		
For this kind of problems, ask for the technical assistance, after having checked the signal by switching off and on the unit		
ERRORS RELATED TO MASTER MICROCONTROLLER		
RAM Memory	Err 32	20
FLASH Memory	Err 34	22
System complemented variables	Err 35	23
Watchdog Timer	Err 33	21
Errors related to power supply section		
+5VDC failure	Err 36	24
HF power supply failure with lower setting value	Err 37	25
HF power supply failure with higher setting value	Err 38	26
HF power supply higher than setting value during activation	Err 97	61
Errors related to serial communication		
Failure regarding serial communication peripheral placed on mother board also code 801463	2 slow sound repetitions; after fast repetitions	3d
Failure regarding serial communication peripheral placed on power supply board also code 801471	3 slow sound repetitions; after fast repetitions	3E
Failure regarding serial communication peripheral placed on double handswitch & nP control board also code 801462	4 slow sound repetitions; after fast repetitions	3F
E ² PROM Communication error	5 slow sound repetitions; after fast repetitions	40
I ² CBUS connection lost (during the use)	6 slow sound repetitions; after fast repetitions	43
Master-Slave connection lost (during selftest phase)	7 slow sound repetitions; after fast repetitions	41

Error description	Error code	Stored code in E ² PROM (select "0d" in CAL mode) ¹
Slave microcontroller failure	8 slow sound repetitions; after fast repetitions	44
ERRORS RELATED SLAVE MICROCONTROLLER		
Failure in a Slave microcontroller internal peripheral	Err 60	3C
Slave microcontroller connection failed	Err 68	44
Errors related to HF power section		
Failure in internal dummy load used for selftest	Err 39	27
Failure in HF output current measurement (lower than expected)	Err 51	33
Failure in HF output current measurement (higher than expected)	Err 52	34
Failure in bipolar output power(lower than expected)	Err 43	2b
Failure in bipolar output power(higher than expected)	Err 42	2A
Failure in monopolar circuit when PURE CUT mode is used	Err 40	28
Failure in monopolar circuit when SPRAY COAG mode is used	Err 41	29
Failed modulation signal for monopolar BLEND-1 mode	Err 44	2C
Failed modulation signal for monopolar FULG mode	Err 45	2d
Failed modulation signal for monopolar SPRAY mode	Err 46	2E
Output peak voltage reading failed (lower than expected)	Err 47	2F
Output peak voltage reading failed (higher than expected)	Err 48	30
Output peak current reading failed (lower than expected)	Err 49	31
Output peak current reading failed (higher than expected)	Err 50	32
HF output power higher than expected	Err 98	62
ERRORS RELATED TO THERMAL PROTECTION		
Power supply thermal protection	Err 27	1b
HF power section thermal protection	Err 28	1C

ENVIRONMENTAL and ATMOSPHERIC CONDITIONS for USE, TRANSPORT and STORAGE

The environmental use and preservation conditions are the following ones:

Temperature (°C): +10 ÷ +40 – Humidity: 30% ÷ 75% - Pressure (hPa): 700 ÷ 1060.

The unit must not be used at less than 30cm from the wall or other objects that can obstruct the ventilation areas, and it must be placed on a trolley or a support (there are two screwed points under the basis for the fixing).

When the unit is not used, it must be kept in a dry place, not dusty. Be aware that no liquid is poured on it.

The environmental transport and storage conditions are the following ones:

Temperature (°C): -40 ÷ +70 - Humidity: 10% ÷ 95% - Pressure (hPa): 500 ÷ 1060.

For the shipment of the unit, we suggest to use the original packaging, or at least a new one that can guarantee the same reliability.

MAINTENANCE and DISPOSAL

The unit must be regularly checked (once per year) by qualified personnel, even better if by the Manufacturer himself. Check always the perfect condition of the accessories, because otherwise their use can be dangerous.

Its disposal must be done according to the different National Rules, by bearing in mind that the accessories (notably the active and neutral ones) get in contact with the patients tissues, and therefore require a special attention.

CLEANING AND STERILIZATION

- Clean the unit with a neutral soap solution, **by taking care that no liquid goes inside**, and then wipe it by a dry cloth. Clean the pedal foot-switches in the same way, or with a cold disinfectant solution (i.e. "Amuchina");

2. Attention: at the moment of the sale the accessories are not sterile.

The following accessories can be sterilized in autoclave (*gloves cycle* at 121°C for 20min / 134°C for 10min), or with cold solutions (i.e. "Amuchina"):

- The electrodes-holder handles and all the active electrodes;
- The cables and the forceps for bi-coagulation, the cables and the forceps (scissors, cannulae) for monopolar coagulation;
- The cables and the monopolar/bipolar electrodes for laparoscopy;
- The handles, the cables and the electrodes for the Argon gas section.

The reusable type neutral electrodes (electrodes and cables) can be sterilized only with cold solutions (i.e. "Amuchina").

During the sterilisation, do not bend too much the connection cables, and wipe perfectly all the parts of the accessories before using them, in order to eliminate all traces of humidity. The best thing to do is to centrifuge them.

TECHNICAL FEATURES

- Electronic generator according to the Standards CEI EN 60601-2-2 (IEC 601-2-2 ed. 1991);
- General mains switch;
- Monopolar and Bipolar working frequency: 435Hz +/- 10%;
- Classification CEI/IEC: Class I - Type CF,
- Classification 93/42 CEI/EEC: IIB;
- Output circuit: "floating", insulated from ground at high and low frequencies, and protected against the use of the defibrillator;
- Mains and Absorption: 230 V ~ 50 Hz - 840 VA, Mains Fuses: T 5 A;
- Monopolar mode with the possibility to connect 1 or 2 electrodes-holder handles (both of them with hand-switches, or just one with hand-switches and the other one with pedal foot-switches);
- Bipolar mode with the possibility to connect 1 or 2 electrodes with pedal foot-switches;
- Bipolar mode in coagulation by pedal foot-switch or (with a particular selection) by automatic start/stop system, according to the conditions of the tissue (impedance sensing); adjustable activation delay from 0.5sec to 5sec;
- Typical working values: 0-30Ω = non activated system; less than 900Ω = start; from 1000 to 1700Ω = stop;
- Functioning memorisation system with 100 memories;
- Pedal foot-switches: 1 or 2;
- Setting: by keys keys for the functioning and the powers (indicated on the display in centesimal scale);
- Functioning control circuit: by twin microprocessor;
- Self-control and self-diagnosis circuit with delivery stop system, acoustic signal, visualization of the fault code;
- Neutral electrode safety circuit with acoustic signal (loud and intermittent), luminous signal (red), and visualization of the code;
- Protection against liquids: standard, non protected enclosure;
- Cooling system: by convection, without fan – Discontinuous functioning: max. 1 hour (10s ON/30s OFF);
- Shelf-life: 5 years;
- Dimensions and weight: (LxDxH) 39x36x17cm – 17Kg;
- Mains cable: 3m long, section 3x1mm².

MONOPOLAR CURRENTS

CURRENT	Maximum Power	Load	Vp/p	Frequency (F) - Crest Factor (CF) - Modulation (M) - Duty Cycle (DT)	Acoustic and Luminous Signals
PURE CUT	350 W	360 Ω	3455 V	F: 435 kHz- CF: 1.6 – M: No – DT: no	Loud tone, Yellow light
BLEND 1	300 W	360 Ω	3595 V	F: 435 kHz- CF: 2.3 – M: 29 kHz – DT: 65 %	Loud tone, Yellow light
BLEND 2	140 W	620 Ω	7600 V	F: 435 kHz- CF: 8.1 – M: 19 kHz – DT: 9 %	Loud tone, Yellow light
ENDO	220 W	360 Ω	1720 V	Alternation 50% of Pure and Blend1 - CF:2.2	Loud tone, Yellow light
FULG	150 W	345 Ω	4700 V	F: 435 kHz- CF: 4.5 – M: 78 kHz – DT: 3.5%	High tone, Blue light
PIN POINT	250 W	250 Ω	3460 V	F: 435 kHz- CF: 2.6 – M: 29 kHz – DT: 50 %	High tone, Blue light
SOFT	280 W	250 Ω	3440 V	F: 435 kHz- CF: 2.5 – M: 29 kHz – DT: 56 %	High tone, Blue light
SPRAY	140 W	630 Ω	7600 V	F: 435 kHz- CF: 8.1 – M: 19 kHz – DT: 9 %	High tone, Blue light

BIPOLAR CURRENTS

CURRENT	Maximum Power	Load	Vp/p	Frequency (F) - Crest Factor (CF) – Modulation (M) - Duty Cycle (DT)	Acoustic and Luminous Signals
PURE	140 W	310 Ω	840 V	F: 435 kHz- CF: 1.5 – M: No – DT: no	Loud tone, Yellow light

BLEND	120 W	310 Ω	990 V	F: 435 kHz- CF: 1.8 – M: 29 kHz – DT: 75 %	Loud tone, Yellow light
MICRO	120 W	100 Ω	440 V	F: 435 kHz- CF: 1.7 – M: No – DT: no	High tone, Blue light
AUTO	See MICRO				
MACRO	120 W	100 Ω	750 V	F: 440 kHz- CF: 1.7 – M: No – DT: no	High tone, Blue light




Conformity EMC/Directive 89/336/CEE: Category A Suggested distances to keep from not vital devices		
Source of the Current RF	Typical Power (W)	Distance (m)
Microcellular telephones CT1,CT2,CT3	0.01	0.4
Mobile telephones DECT, Wireless devices (modems, LANs)	0.25	2
Mobile telephones (USA)	0.6	3
Hand mobile telephones (GSM, NMT, Europe)	2	6
(DECS 1800)	8	11
Walkie-talkie (police, firemen , protection, maintenance)	5	9
Bag mobile telephones	16	16
Mobile radio (police, firemen, protection)	100	40
For broadcasting stations which use frequencies less than 800MHz, the distance can be established by using the equation: A: $d = 4\sqrt{P}$ For broadcasting stations which use frequencies between 800MHz and 2.5GHz, the distance can be established by using the equation: B : $d = 2.3\sqrt{P}$ P = Nominal power of the transmitter in watt (W), established by the manufacturer.		

MONOPOLAR CURRENTS




(Variation of the output power from 50 to 2000Ω by selecting 100% and 50% of maximum power)

(Increase in the output power with nominal load according to the power adjustment) (Increase of the Voltage according to the power adjustment).




Corrente / Current / Courant / Corrente "PURE"

Corrente / Current / Courant / Corrente "BLEND 1"




Corrente / Current / Courant / Corrente "BLEND 2"




The better way to operate TM




Corrente / Current / Courant / Corrente "ENDO"

Corrente / Current / Courant / Corrente "FULG"




Corrente / Current / Courant / Corrente "PINPOINT"




The better way to operate TM

Corrente / Current / Courant / Corriente "SOFT"

Corrente / Current / Courant / Corriente "SPRAY"




BIPOLAR CURRENTS




(Variation of the output power from 10 to 1000Ω by selecting 100% and 50% of maximum power)

(Increase in the output power with nominal load according to the power adjustment) (Increase of the Voltage according to the power adjustment).




Corrente / Current / Courant / Corrente "PURE"




Corrente / Current / Courant / Corrente "BLEND"




Corrente / Current / Courant / Corrente "MICRO"

Corrente / Current / Courant / Corriente "MICRO AUTO"

Corrente / Current / Courant / Corriente "MACRO"


ACCESSORIES

ELECTRODES-HOLDER HANDLES (nominal voltage: 4000 V)

Electrodes-holder handle (electrodes with stem Ø 2,3/2,4 mm) with cable 3.5m long;

Electrodes-holder handle (electrodes with stem Ø 2,3/2,4 mm) with hand-switches, reusable 100 times.


ACTIVE ELECTRODES STANDARD SERIES (see next page for details)

(Nominal voltage: 4000 V)

- E 1 - Straight knife electrode
- E 5 - Straight needle electrode, thick type (2 pcs.)
- E 7 - Straight needle electrode, thin type
- E12 - Straight ball electrode Ø 2.5 mm
- E14 - Straight ball electrode Ø 4 mm (2 pcs.)

NEUTRAL ELECTRODES

Stainless steel neutral electrode for adults (bodies with weight higher than 15 Kg), cable 2.5m long;

Flexible rubber neutral electrode for adults (bodies with weight higher than 15 Kg), cable 2.5m long;

Flexible rubber neutral electrode for children (bodies with weight from 5 to 15 Kg), cable 2.5m long.

PEDAL FOOT-SWITCHES

94-933: Twin pedal foot-switch;

34-9331: Twin pedal foot-switch for bipolar mode only.